




Data Analysis Software DS-5 (Ver. 2.10)



KONICA MINOLTA SENSING, INC.

INSTRUCTIONS FOR USE

DS-5 may be used with
KONICA MINOLTA Pulse Oximeters:
Pulsox-3i/ 3iA/ 3Si/ 3Li and Pulsox-300i

Contents

FOR YOUR SAFETY AND THAT OF YOUR PATIENTS	4
Information on safe use	4
Operation	5
SETTING UP A PULSOX-3 SERIES OR PULSOX-300I DEVICE PRIOR TO, & MAKING A RECORDING	5
Starting the DS-5 Programme	5
'Download' a recording from a Pulsox 300-i oximeter	6
Setup: Date and Time	7
Setup: Body Number	8
Setup: Duration of the recording	8
Data File Operation: Clear Memory	8
Make the Recording	9
On return, after the recording has been made	9
Data File Operation: Download	9
Viewing previously downloaded data	12
Viewing the 'Chart' screen	13
Navigation through the recording	14
Showing the SpO2 dips and Pulse Rate rises	14
Changing the Compression 'Method'	14
Revert	14
ANALYSIS: RUN AN ANALYSIS AND MAKE A REPORT	15
Analysis Customisation	15
Report Setup	15
SpO2 Analysis Setup	16
Pulse Rate Analysis Setup	17
Graph	19
ANALYSES:	20
Analysing data	20
The Analysis screen	21
Ending the programme	21
ADVANCED FEATURES	22
Patient Details: Entering and Editing	22
Editing the Labels	25
MONTAGES	26
MARKERS	27
EVENTS	28
PRINTING	30
Print Window	30
Print	30
Print Setup	30
The Print Setup dialogue box is supplied by the installed printer.	30
SAVING, COMBINING, IMPORTING AND EXPORTING DATA	31
Open and Close	31
Close Workpad File	31
Save	31
Save As	31
Export	31
Import	33
Combine Files	34
Analyse and Dump (Export) Patient information	34
Combine Patient Information Files	35
Group Auto Operation	35
Delete (a recording, analysis) and Delete Workpad	36
Properties: (Information about the recording)	36
Window Menu	36
WINDOW CHOICES	37
STANDARDS	38
Safety and effectiveness	38
Indications for use	38
User Responsibility	38
Warnings and contra-indications	38
Standards compliance	38

Symbols 38

Specifications..... 38

Manufacturer's Address..... 38

APPENDIX A: DUMP DATA..... 39

APPENDIX B: KONICA MINOLTA PULSOX-3 SERIES OXIMETERS: 3I/ 3IA /3Si/ 3LI 40

DS-5 Software v2 © 2007 KONICA MINOLTA SENSING, INC.

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the software requires full understanding and strict observation of these instructions.
The software is only to be used for the purposes specified here.

Maintenance

The software installation must be carried out and serviced by trained service personnel. Reinstallation and general maintenance may only be carried out by trained service personnel.

Accessories

Only the devices listed in the Devices List may be used with this software.

Safe Connection with other electrical equipment

Electrical connections to equipment which is not listed in these 'Instructions for Use' should not be made.

Liability for proper function or damage

The Warranty/ Software license agreement is included on the delivery CD.

KONICA MINOLTA SENSING, INC. does not warrant that the software will be uninterrupted or error free.

Information on safe use

The software should only be used by trained users of Microsoft Windows.

DS-5 may only be used in the home or hospital.

DS-5 is intended for use with CE marked devices

DS-5 is designed to view, analyse, report and print data recorded by these devices. Inherent in its design is the capability to review or fully disclose the 'raw data' obtained.

DS-5 can only deal with the stored results of the measurement by the device. It does not influence the measurement, it is reviewed after the fact and has no direct relevance to patient safety.

OPERATION

Windows 7

Windows 7 (and Windows Vista) has replaced:

C:\Documents and Settings\All Users\Application Data\

folder with:

C:\ProgramData\

If you are not able to see the C:\ProgramData\ folder go to folder properties and enable show hidden files and folders.

Setting up a Pulsox-3 series or Pulsox-300i device prior to, & making a recording

Follow the instructions for use in the Device Manual.


In particular, prior to making a recording set date and time, recording requirements and delete any existing recordings as indicated below.

Starting the DS-5 Programme


- Choose the DS-5 icon on the desktop or DS-5 from the Start | Programs menu



Initial 'Splash' screen as programme opens



DS-5 showing shortcut icon to Download the Pulsox-3 series devices




DS-5 showing shortcut icon to Download the Pulsox-300i device

‘Download’ a recording from a Pulsox 300-i oximeter


(See ‘APPENDIX B: KONICA MINOLTA PULSOX-3 SERIES OXIMETERS: 3i/ 3iA/ 3Si/ 3Li’ to download a recording from the Pulsox-3i series.)

- Click on ‘Communication menu’ or click on the Download icon appropriate to the device in the ‘Toolbar’. As the mouse passes over the icon, the device name will appear in the ‘yellow box’



Download Pulsox-300i dialogue box before the Pulsox-300i is connected/ or turned on

The recorder is connected to the computer using the USB interface. It is necessary to set up this interface when first used, for which see the device User Manual.




Download Pulsox-300i dialogue box when the Pulsox-300i is connected and turned on

The dialogue box with items in black shows that communication with a recorder (Body No. 42) has been established.

- *Setup* and *Data File Operation* allow setting of the recorder parameters and data etc.

Setup: Date and Time

- Click on the Set Clock button in the Download Pulsox-300i dialogue box
- The Set clock information box shows the data and time obtained from the pc.




Set Clock information box

- Click on the Set button to send the information to the recorder.

Setup: Body Number

- Click on the Set Body No. button in the Download Pulsox-300i dialogue box




Set Body No. dialogue box

- After the device has been interrogated, a dialogue box will show the Body No., which can be changed and Set

Setup: Duration of the recording

- Click on the Set Auto Power Off Time button in the Download Pulsox-300i dialogue box



Set Auto Power Off Time dialogue box

- You may choose the duration of the recording from switch-on. To disable this feature tick 'Do not set auto power off'

Data File Operation: Clear Memory

NOTE: It is normal to start a new recording with the memory cleared (deleted).

Click on the Clear Memory button in the Download Pulsox-300i dialogue box

A warning box appears before the memory is cleared

**Clear Memory warning box**

- Clicking on Yes will delete all data from the recorder

Make the Recording

Normally at this point the recorder is removed and attached to the patient. For this, follow the instructions in the Device Manual.

On return, after the recording has been made

When the device is returned from the patient you are ready to copy or 'download' the recorded data into the computer.

Restart the program.


Click on the 'Communication menu' or click on the Download icon appropriate to the device in the 'Toolbar'.

Connect the recorder using the USB interface.

The dialogue box with items in black, shows that communication with a recorder has been established.
(as shown in **'Download Pulsox-300i dialogue box when the Pulsox-300i is turned on'**)

Data File Operation: Download

To download the data first click on the Download button.



Download Pulsox 300i Series dialogue box with several recordings on the recorder, file #2 chosen and an automatically generated filename

The Download Pulsox 300i series dialogue box will then show the recordings on the recorder.

To download recordings from the recorder into the computer

- Highlight the file using the left mouse button
- Click on the Download button

Show All Files/ Hide Short Files Button

Short recordings of less than 300 seconds duration may be hidden from view or shown with this button.

File Selection

Highlight one or more recordings. You can select two or more recordings by pressing the <Shift> or <Ctrl> key simultaneously with the left mouse button.

Note: If you select two or more recordings DS-5 will automatically merge the recordings into one file with zero values added to fill the gaps between the recordings. These files may be very large (up to 14 days) so if you do not want to merge recordings they should be downloaded one at a time.

Location of downloaded file / Browse

'Location of downloaded file' shows where the selected recording will be saved on the computer. By default this is the program data directory.

To change where the file will be saved click the 'Browse' button and select a folder or create a new folder.

File Name

A unique File Name is automatically generated and the default location of the file is indicated.
A different filename may be entered and the location may be changed with the Browse button.

Surname / ID Number

You may enter the patient surname/ID to be saved in the file. If you leave the patient name field empty it will automatically be saved with a blank patient name, which could lead to problems in locating the file in the future.

Details

This gives access to a dialogue box to allow patient details to be entered. See 'Patient Details' section for the details.

Download Channels

Choose the channels to be downloaded by marking the check boxes.

Auto print report after downloading

You may select 'auto print report after downloading', to download, analyse and print the report in one operation.

Auto create PDF report after downloading

You may select 'auto create PDF report after downloading' to download, analyse and create the PDF report in one operation.

Auto create patient information file after downloading

You may select 'Auto create patient information file after downloading' to download, analyse and create the patient information file in one operation.

Report Setup


Before downloading you may use the Analysis setup dialogue box to define the parameters used in the analysis. See 'Analysis: Run an Analysis and make a Report' section for the details.


Downloaded

When the Download button is clicked the highlighted recordings will be transferred to the computer. The data is saved by this action but it is not deleted from the recorder. The number of times a file has been downloaded is shown.

After the data has been transferred, the first few minutes will be shown on the 'Chart Screen'.





Viewing previously downloaded data


- Click on the file open icon 
- The Open dialogue box shows the recordings which are stored in the data directory.
TIP: If it does not show the 'Data' directory in the Current Open Directory use this drop-down menu to navigate around the computer.
- Highlight the required recording: the yellow box will show some identification data that has been entered about the subject - and click on Open or double click on the patient name to open it to the Chart Screen.



The Open dialogue box

There may be a **cross**, a **filled circle**, or a **dash** in front of the Patient names in the **Files** box.

-  indicates that there is no saved analysis file for the recording.
-  indicates that an analysis of this recording or markers has(ve) been made and saved as a Workpad file.
-  appears in front of a subject with a Workpad analysis(es) which follows, preceded by a 




Double clicking on a cross will display saved Analysis or 'Workpad' files for this patient and the cross becomes a dash.

Then select the required analysis workpad.

NOTE : Opening a workpad file will also open the corresponding recording.

Viewing the 'Chart' screen

When a recording is first opened, whether after download or in Review, the first interval will be shown on the screen.



File | Open | Sample selected


The Chart screen consists of

- **Title Bar** - shows the programme name, patient name and filename.
- **Menu Bar** - has a series of functions grouped in menus.
- **Tool Bar** - has chosen items from the Menu Bar.
- **Ruler** - Marks between a pair of moveable cursors the interval shown in the Chart Data Window, allowing you to see how it relates to the recording.
- **Chart Window** - showing the channels of the recording, SpO2, Pulse Rate, etc with the values at the cursor time appearing at the right hand side of each channel. At the top of the trace channel are shown the alarms and events generated by the device. Placing the cursor pointer over the marking will reveal the alarm/ event type
- **Status Bar** showing the interval duration, real date and time at the cursor, elapsed time, value at the point of the arrow cursor, compression rate and montage if applied.



Change the appearance of the Chart Window

Move the mouse over the graph you wish to modify and click the **Right mouse button**. The menu functions:

- **Amplitude** -change the vertical axis of the trace over which the mouse is placed using mouse selected, movable cursors, eg Saturation shown here. After completing the axis modification action, close the box with a click on the top left hand square
- **Exclude** excludes the trace
- **Exclusive** shows the trace **exclusively**.
- **Parameter** allows numerical manipulation
- **Fit To Scale** uses the raw data rather than the data as viewed in the current **Method**.





Navigation through the recording


Move through the recording page by page using the **Page Right** (or **Page Left**) icons  in the Tool Bar or by clicking on the right (to go forward) or left (to go back) on the trace. The Pan icons  step through the trace at the page display duration set in menu item View | Page | Pan Rate.

Display Mode

The display mode can be changed by several means.

- For a different Interval (time scale) click on the Interval icon  in the Tool Bar.
- To see the whole recording click on the Full Study  in the Tool Bar
- For a different amplitude use the mouse right hand button, highlight and choose Amplitude to change the vertical axis by moving the cursor(s) - or choose Parameter to set the values numerically. (Amp. is the value at the top of the axis, Offset is that at the bottom)
- The user may wish to take advantage of the ability to create 'Montages': combinations of screen position, amplitude range, colours, etc) for favourite repeated use. Choose **M** on the toolbar to use one of the 5 default montages. (These montages can be changed and saved using the menu item View | Edit Montages)
- To show a grid select View | Show Grid (change colour of lines with Edit | Gridlines colour).

Showing the SpO2 dips and Pulse Rate rises

Select the Toolbar icon  to show the dips and rises as previously defined (see the Analysis section).

Changing the Compression 'Method'

It is not possible to display all the recorded data points if a long time interval is displayed on the screen. One point on the screen (called a pixel), will represent a 'package' or group of data. The 'pack rate' shown in the Status Bar indicates the number of data points represented by one pixel. The different ways to represent this package on the screen are:

- **Sample** : displays the first value of each package
- **Mean** : displays the average of all the values in the package
- **Minimum** : displays the lowest value in the package
- **Maximum** : displays the highest value in the package
- **MinMax** : displays the lowest and the highest values in the package.
- **Envelope** : displays a line joining the lowest and the highest values in the package

The **Method** can be changed for an individual channel with the right mouse menu item Parameter.

Revert

Go back to the way the recording appeared originally, by clicking on the Toolbar icon **M** and then **O: Original montage at last save**

Analysis: Run an Analysis and make a Report


If the recording is acceptable throughout, the next task is to analyse the data and print a report. However, there may be periods that you wish to exclude, for which see the Markers section below.

Analysis

Analyses are selected in the Analyse menu, or with certain icons, eg. to make a comprehensive analysis,

select **Analysis Report** directly with the toolbar icon 


Analysis Report: The Analysis Report consists of several pages, chosen by ticking in the Analyse | Analysis Setup | Report | Report Contents

To create a Adobe Acrobat .pdf version click on Analysis Report to pdf 

Or: Select one or more of the analysis options from the **Analyse** menu and its sub-menus.

- **Analysis Setup** to define the parameters used in the analysis (see below)
- **Analysis Report** is a full report as defined in the Analysis Setup
- **Analysis Report to PDF** is a full report in PDF format
- **Analyse** options produce individual analysis pages appropriate for sleep study, day, or walk test
 - Analyse SpO2 (Sleep study)**
 - Analyse SpO2 (Sleep-80 dips)**
 - Analyse SpO2 (Daytime)**
 - Analyse Pulse**
 - Walk Test**
- **Analysis Summary ()** produces a summary page of saturation and pulse rate analysis.
- **Plot selected channels** show graphically 1 hour per line, 12 hours per page (default).
- **Plot selected channels (12 or 24 hrs)** are 2 graph per page plots
- **Auto print report** analyses and prints a set of plots, analyses and summary pages as defined in Analysis Setup.
- **Copy Analysis Summary () to clipboard** to paste the summary into Word etc
- **Display analysis using edited events** to reanalyse after creating or modifying events

Analysis Customisation

Choose **Analyse | Analysis Setup**, or the **Analysis Setup** toolbar icon 

Report Setup

You may select pages and options for the Analysis Report

Report Contents:

Select the pages to be included in the report:

Report Sub-Heading

Enter the hospital or clinic details, telephone etc. for printing on the Report
Line 1 is printed in a larger font.

Tick '**Show extra patient details on report**' to show on the Patient Details the values entered alongside the defined labels in the Extra Patient Details sheet.

Report Options

- **Plot graphs using channel colours** Select by ticking (un-tick for black & white printers)
- **Hours to plot per page.** Selected channels are those selected in [Analysis Setup] [Graph] 'Select Plot Channels' window.

Report Markers

Select **Full Study**, **Within Markers** or **Outside Markers**. **Within Markers** and **Outside Markers** are effective only if user insert markers.

Analyse | Analysis Setup | Report

SpO2 Analysis Setup

Dip Definitions

- **SpO2 cutoff for analysis**
Select the cutoff value: the level below which the data will not be included in the analysis.
- **Dip definition (1) / (2) / (3)**
Enter the saturation dip (in %) value: greater than which the oxygen saturation needs to drop from the last peak value to constitute a dip. To enable the calculation of definitions 2 and 3 tick the boxes left of Dip definition (2) and/or Dip definition (3).

Frequency Graph

- **Maximum time to plot**
Select the value used as the maximum value for the vertical axis of the frequency histogram plot.
- **Auto adjust frequency graph**

Set the vertical axis to automatically adjust to the maximum rather than using the Maximum time value

Chart Screen

- **Plot SpO2 mean**
Tick this option to plot the mean value as a dotted horizontal line on the chart screen. (the mean value calculation excludes values below the Cutoff point).

Other

- **Time spent SpO2 < xx%**
Tick this option to calculate the time under a required value
- **Moving Average Time**
Set the period over which the SpO2 is averaged Moving average can be applied for Pulsox-300i only. Moving average is not applied to calculating sleep apnoea index.

Select Display

- **ODI** Displays ODI (Oxygen Desaturation Index) results on analysis plots
- **BMI** Displays BMI (Body Mass Index) results on analysis plots

- **Delta Index** Displays Delta Index on SpO2 analysis (Sleep Study) report

Dip Counting Method

Set the times used in the analysis

Results of applying new parameters

- **Dips:** the total number of desaturations defined in Dip definition (1)
- **Dips/hr:** total number of desaturations divided by the analysis period in hours
- **Mean Nadir:** the mean of the lowest values found within each dip
- **Mode:** Full Study, Within or Outside Markers

Apply: Test and display the analysis immediately (these results are only visible, not saved)

Reload: For reinstating theSpO2 analysis setup defaults.

Analysis Setup

Report | SpO2 | Pulse Rate | Graph

Dip Definitions

SpO2 Cutoff for analysis: 25

Dip definition (1): 4

☒ Dip definition (2): 3

☒ Dip definition (3): 2

Time spent SpO2 < 88 %

Select Display:

☐ ODI ☒ BMI ☐ Delta Index

Frequency Graph

Maximum time to plot: 2:00:00

☒ Auto adjust frequency graph

Chart Screen

☒ Plot SpO2 mean

Moving Average Time (Pulsox - 300i only)
(not applied to calculating sleep apnoea index)

☒ 3 (sec) ☐ 6 (sec) ☐ 12 (sec) ☐ 24 (sec)

Dip counting method

SpO2 drop >= a

SpO2 rise >= b

c1 <= actual time for drop <= c2

actual time for rise < d

a : dip definition

b : rise after drop

c1 : Minimum time to drop in seconds

c2 : Maximum time to drop in seconds

d : Maximum time to rise in seconds

b = a

4

120

20

Results of applying new parameters

	>=4%	>=3%	>=2%
Dips:	30	64	153
Dips/Hr	1.25	2.67	6.38
Mean nadir:	89.47	90.45	91.86

Mode: Full Study

Reload Apply

Save Cancel

Analyse | Analysis Setup | SpO2

Pulse Rate Analysis Setup

Rise Definitions

- **Pulse cutoff for analysis**
Select the cutoff value: the level below which the data will not be included in the analysis.
- **Rise definition**
Enter the amount (in bpm) by which pulse rate needs to increase to constitute a rise. Three possible rises can be calculated. To enable the calculation of definitions 2 and 3 tick the boxes left of Rise definition (2) and/or Rise definition (3). The default Rise definition (1) is 6 bpm.

Plot Options

- **Maximum time to plot**
Select the value used as the maximum value for the vertical axis of the pulse rate frequency graph.
- **Auto adjust frequency graph**
Set the vertical axis to automatically adjust to the maximum rather than using the Maximum time value
- **Plot Pulse mean**

Plot on the chart screen the mean value (the mean value calculation excludes values below the **Pulse cutoff for analysis**).

The screenshot shows the 'Analysis Setup' dialog box with the 'Pulse Rate' tab selected. The 'Rise Definitions' section on the left includes a 'Pulse cutoff for analysis' set to 25, and three checked 'Rise definition' options (1, 2, and 3) with values 6, 7, and 8 respectively. The 'Frequency Graph' section has a 'Maximum time to plot' of 02:00:00 and a checked 'Auto adjust frequency graph' option. The 'Chart Screen' section has a checked 'Plot Pulse Mean' option. On the right, the 'Results of applying new parameters' section is empty, with labels for 'Rises:', 'Rises/Hr:', and 'Mode'. At the bottom right are 'Reload' and 'Apply' buttons, and at the very bottom are 'Save' and 'Cancel' buttons.

Analyse | Analysis Setup | Pulse Rate

Results of applying new parameters

Apply: Test and display the analysis immediately (these results are only visible, not saved)


Reload: For reinstating the Pulse Rate analysis setup defaults.

- **Rises:** Number of pulse rises defined in Rise definition (1)
- **Rises/ Hr:** Number of pulse rises divided by the time in hours
- **Mode:** Full Study, Within or Outside Markers

Graph

Graph contents

Select the channels to be plotted and their lowest and highest values



Analyse | Analysis Setup | Graph

Analyses:

There are two plots shown on the analysis pages. On the left is the Interval Frequency Histogram (the vertical axis is set in Setup, either as a fixed value e.g. 2.00.00 (hh.mm.ss), or as auto-adjust.) On the right hand side is a Cumulative distribution graph.

The horizontal axis of the saturation graph and the horizontal axis of the Pulse Rate are set in the Analysis Setup | Graph dialogue box.

Analysing data

General parameters for saturation, and pulse rate analysis

- **Mean** (the average) of the analysed values
- **Median** the value below which 50% of the points fall
- **Duration** Duration of the Full study
- **Marker Mode** either **Full Study**, **Within** or **Outside** (the Markers)
- **Analysis Interval** **Duration of the full study, within or outside markers**
- **Total time analysed for ‘ ‘** the time interval which has been analysed for ‘ ‘

Parameters for SpO2 Analysis

- **5th/ 95th percentile** the lowest and highest 5% of the data (also drawn on the plots)
- **Time spent SpO2 < 95% etc** the period (time interval) and % of analysed time
- **Time spent SpO2 < xx%** as defined in 'Analysis Setup'

Parameters for Sleep Apnoea Index

- **Dips** dips calculated from a moving baseline, the % required to constitute a dip may be defined by the user (see above).
- **Dips /Hr** number of dips/ analysed period in hours
- **Mean nadir** mean of the lowest saturation level in dips.
- **Min Dip SpO2** lowest saturation level in a dip.
- **Delta Index** a measure of variability in the data: see below

Delta Index

This is an index of variability in the oximetry trace. In the references, oxygen saturation values stored at 12 second intervals from oximeters were summed, and the absolute sum, divided by the total number of samples was calculated as the **Delta Index**.

The oximeters used by the authors store the lowest value in the preceeding 12 second period unless it is lower than 80% of the mean of the rest of the values. ('Lowest value rule')

For other oximeters, the data is mapped to a 12 second rate using digital interpolation, with the 'lowest value rule' applied. To the data points: a 'lead in' time of 12 seconds is used.

A second, oximeter specific 'delta index' using the sampling rate of the oximeter e.g. 5 seconds, is also calculated.

References: Levy P., Pepin JL., Deshaux-Blanc C., Paramelle B., and Brambilla C, 'Accuracy of Oximetry for Detection of Respiratory Disturbances in Sleep Apnea Syndrome' Chest 1996;109: 395-99.
Pepin JL., Levy P., Lepaulle B, Brambilla C and Guilleminault C. 'Does Oximetry Contribute to the detection of Apneic Events? Mathematical Processing of the SaO2 Signal.' Chest 1991 99 1151-57.


Parameters for the Pulse Rate Analysis

- **Mode** the most frequently found pulse rate
- **Standard Deviation** (the square root of the variance) -a variability index


Parameters for Pulse Rises

- **Pulse Rises >xx** pulse rises above the background.
- **Rises/ hr** the average hourly pulse rise rate
- **Max Pulse Peak** highest pulse rate level in a pulse rise

The Analysis screen

After running an analysis the first page of the report is visible. From here all pages can be displayed, by clicking the toolbar icon **All Pages**  or by the right hand mouse button.

The left mouse button can be used to select and zoom in on a particular (mouse selected) page, and the right mouse button to return.



Example 'All Selected Pages' of Report

Ending the programme

To exit the programme click File | Exit

This closes the opened documents. A 'Save Changes?' box will appear for each changed, unsaved document, and the application will then be closed.

Advanced Features

Patient Details: Entering and Editing

To change the patient name stored with the current recording, (as displayed on the top of the Chart screen) choose **Edit | Patient Details** which gives a four tab dialogue box, the first tab showing Patient Details

Edit | Patient Details | Patient Details

- Up to 40 characters are available for the Surname.
- A predefined comment can be entered using Comment Selection.
- Plain text may be entered in the Comment Input box whenever the dialogue box is opened.
- Body Mass Index is calculated from the metric Height and Weight data
- A subset of this data appears in the 'Details' (yellow box) in the 'Open' box

Fixed Form Click on the Fixed Form button to amend the ten available predefined fixed form comments

Save using the **Save** button.

The surname is the name displayed in the title bar. It is updated in the Open dialogue box after being saved.

Walk Test

Using the second tab, enter the Borg Scale for dyspnoea and fatigue measured at the start and end of the Walk Test, the supplemental oxygen used and distance walked

The screenshot shows a software window titled "Edit patient details" with four tabs: "Patient Details", "Walk Test", "Medical Details", and "Extra Patient Details". The "Walk Test" tab is selected. The form contains the following fields:

- Borg Scale**
 - Baseline**
 - Dyspnea: Dropdown menu showing "1: Very slight".
 - Fatigue: Dropdown menu showing "2: Slight (light)".
 - End of Test**
 - Dyspnea: Dropdown menu showing "5: Severe (Heavy)".
 - Fatigue: Dropdown menu showing "7: Very severe".
- Supplemental oxygen flow during the test**: Text input field showing "1.0" with the unit "(L/min)".
- Total walking distance**: Text input field showing "234" with the unit "(m)".

At the bottom right of the window are "Save" and "Cancel" buttons.

Edit | Patient Details | **Walk Test**

Medical Details

Using the third tab enter **Medical Details**, ie **Referring Physician, Health Authority, GP details**

Edit patient details

Patient Details | Walk Test | **Medical Details** | Extra Patient Details

Referring Physician: Dr Fred Evans Health Authority: S Oxon

GP Name: Simon Young

GP Address1: The Surgery

GP Address2: Bletchington Road

GP City: Oxford

GP County: Oxon

GP Post Code / ZIP: OX5 8TY

GP Tel: 01993 456789

Number of times changed: 6 Date last changed: 2007/04/20

Fixed Form

Save Cancel

Edit | Patient Details | Medical Details

Extra Patient Details

The fourth tab: **Extra Patient Details** has 12 boxes, with configurable text lines, set in Edit | Label details. Here the configurable text are Adverse Effects, Sleep disruption etc

Edit patient details

Patient Details | Walk Test | Medical Details | **Extra Patient Details**

Custom Label 1 data 1

Custom Label 2 data 2

Custom Label 3 data 3

12 empty text boxes for configurable text

Edit Label

Save Cancel

Edit | Patient Details | Extra Patient Details

Editing the Labels

In the Edit | Patient Details | Extra Patient Details ‘Edit Label’ button enter your own labels. The new labels may be applied to the current recording (file) and become the default if changes are saved on exit from the recording. It is also accessible by Edit | Label Details.

Label Details

Fixed Labels

Surname:

Address Line1:

Health Authority:

Comments:

Firstname:

Address Line2:

Referring Physician:

Title:

City:

GP Name:

Date Of Birth:

County:

GP Address Line1:

Age:

Post Code / ZIP:

GP Address Line2:

Height:

Tel. Home/Work:

GP City:

Weight:

GP County:

Sex:

GP Post Code / ZIP:

ID Number:

GP Telephone:

Necksize:

Epworth score:

Definable Labels

New Label	Field	New Label	Field
Custom Label 1	(20 chars)		(20 chars)
Custom Label 2	(20 chars)		(20 chars)
Custom Label 3	(20 chars)		(20 chars)
	(20 chars)		(20 chars)
	(20 chars)		
	(20 chars)		
	(20 chars)		
	(20 chars)		

OK

Apply To File

Cancel

Edit | Label Details

Montages


The Chart screen can be changed to have traces with different combinations of screen position, amplitude range, colours, etc. by choosing a **Montage**. There can be five montages defined, M1 to M5.

To apply a montage

- Click the montage icon **M** on the Tool Bar, highlight and click the required montage

Select **Edit | Edit Montages** to edit a montage

The Montage tab shows the montage (1 - 5) currently being edited. The current definitions are retrieved and displayed in the dialogue box. The right hand side values (Heading etc) are those of the highlighted channel.



Edit | Edit Montages


- Clicking the '**Current**' button loads the settings from the currently opened screen into the selected montage.
- Clicking **Save** causes the montage to be validated and then saved.
- **Cancel** abandons any changes and closes the dialogue box.

Note The Montage data is saved to Montage.ini.

Markers


Select **Marker | Insert All Channels Marker** or its toolbar icon to mark with the mouse cursor (click, drag and drop) periods for inclusion or exclusion from saving, export, and analysis.

Select **Marker | Insert Single Channel Markers** or its toolbar icon to mark periods on a single channel



Marker |

Select **Marker | Markers...** to numerically enter new All Channel markers, e.g. from a patient diary, and existing ones to be edited and annotated. Markers previously entered graphically can also be annotated and edited.



Marker | Markers | New Marker

Select **Marker | Goto (Previous/ Next)** options when there are markers already inserted, use the **Marker | Show Markers** to show or hide markers with **Marker | Show annotations** for annotations.

Select the right hand mouse button when in normal viewing mode (no Marker insert selected) to numerically edit or delete both All Channel and Single Channel markers.

Events

On the **Ruler**, after an analysis has been performed, the events which have been detected can be seen.






Chart screen Ruler

In the example above the desaturations are shown in red, and the pulse rate rises (which were analysed later) in blue below the desaturations on the ruler.


On the Chart screen select **Event | Insert/Edit Events** or the toolbar icon  to mark or edit analysis events. Events appear after analyses shown as horizontal bars below the appropriate channel and can be added or modified by the user:

Select **Event | Display/ edit event definitions** to view/ edit/ add event definitions.




Event | Display/ Edit Event Definitions

To create a new event type, enter an event name and click on the Add button. This new event will appear in the 'Events' box, and it should then be highlighted and associated with one of the 'Available Channels' by crossing the appropriate box. The Event Colour may be chosen. It can be included or removed from the default set by clicking on the Delete button(s).



Events marked under the channel on the Chart screen

Click on the Insert/ Edit events icon  to show the events as bars over the Chart which can be edited using the mouse.





Chart after the Insert/ Edit events icon has been clicked

The event information (yellow box) appears when the mouse is positioned over the bottom bar of an event and held steady. The events can be moved, or their boundaries altered using the left mouse button, dragging and dropping, or deleted using the right mouse button on the main body of the event.



Insert Event palette, Dip highlighted for inserting into chart

New events can be inserted in an area where there are no overlapping events by highlighting an event type on the Insert Event palette, and marking them on the Chart with the left mouse button.

When events have been changed, the data can be reanalysed with **Analyse | Display analysis using edited events.**

Printing

Print Window


Select **File | Print Window** to print the current page

or Click on the **Print Window** icon  in the toolbar

The options in Print Setup are used and it is recommended that you check these first.

Print

Select **File | Print** to print the current active document




File | Print

Select the options

- In the **Views** box the channels to be printed by clicking in the box (a cross means that the signal will be printed).
- In the **Analyses** box the analysis(es) to be printed.
- In the **Intervals** box the time interval which will be displayed on each page
- **Full study** to print the full study,
- **Display Interval** to print the interval displayed
- **Time per Page**, use the scroll bar, and highlight the interval to be printed on each page.
- The **Printer** button accesses **File | Print Setup**

Print Setup

The Print Setup dialogue box is supplied by the installed printer.



File | Print Setup

Saving, Combining, Importing and Exporting data

Most of the File menu items become available when a recording is open, or has been modified, but are otherwise greyed out.

Open and Close

Opening has been covered above.

To close the document click on the option **File | Close**,

NOTE: Close all the opened documents, with **Window | Close All**.

Close Workpad File


The Workpad is a document, which holds the analyses, patient details, events etc. It is associated with a recording, but is saved separately. Workpad files have extensions wp0, wp1 etc.

Save

To **Save** after a change has occurred click on **File | Save**

Save As

Use **Save As** to save all / part of a currently open file as another file.



File | Save As

Title Enter a new patient name if required

New File delete the asterisk * before .dat and insert the new filename

If you select a filename of a file that already exists, you will be prompted to overwrite the existing file with the current file.

Directories - shows the directory where you will store the new file.

The box allows navigation.

Channels to save check one or more channels

Markers

If there are no markers, **Within markers** and **Outside markers** will be greyed.


If markers are present, select to save the data for the **Full study**, **Within markers**, or **Outside markers**.

If you select **Within markers** and there is more than one marker present the data between markers will be padded with zeros

If you select **Outside markers**, the data within the markers will be padded with zeros.

Export

To 'export' data, or a subset of the data in a format that can be used in a spreadsheet, or another analysis package.



File | Export

List Files of Type:

- **ASCII or 'text format'** arranges the data in columns according to channel, separated with a single space. The first part of the file forms a header, see below, the first few lines (depending on how many channels were saved) contain information about each channel, chan1 is the data in column one, chan2 column two and so on.

Chan 1: Saturation Sampled at every xx seconds
 Chan 2: Pulse Sampled at every xx seconds
 97.00 86.00
 97.00 85.00 A
 97.00 86.00 etc.

Following the header is the actual data, the number of columns will vary depending on how many channels were saved. Alarm mark is attached to each data.

Alarm mark	Explanation
C	Probe not connected
I	Light improper (data precision not guaranteed)
L	Light improper (impossible to measure)
p	Low perfusion (data precision not guaranteed)
P	Low perfusion (impossible to measure)
H	High perfusion (impossible to measure)
A	Artefact
U	SpO2 range under (data precision not guaranteed)
u	Pulse rate range under (data precision not guaranteed)
o	Pulse rate range over (data precision not guaranteed)

- **DIF (Data Interchange Format)**, commonly used in spreadsheet programs, such as Microsoft Excel™.
- **European Data Format**, usually has a .rec extension and was defined by Kemp et al 1992 'A simple format for exchange of digitised polygraphic recordings E & Clin Neurophysiol 82 (1992) 391-393. http://www.hsr.nl/edf/edf_spec.htm
- **Directories** - shows the directory where you will store the new file.
- The box allows navigation.
- **Channels** check one or more channels
- **Markers**

If there are no markers, **Within markers** and **Outside markers** will be greyed.

If markers are present, select to save the data for the **Full study**, **Within markers**, or **Outside markers**.

If you select **Within markers** and there is more than one marker present the data between markers will be padded with zeros

If you select **Outside markers**, the data within the markers will be padded with zeros.

- **Continuity** - normally each data record represents a discrete time interval and time can be calculated using a simple count of the records, however, if you export data from within markers or outside markers there will be gaps in the data. Continuity deals with these gaps; the **Keep time** option means that the gaps are filled with zeros and thus the time interval is kept, the **Lose time** option means that the gaps are not filled and thus the time interval is lost.

Import

IT IS RECOMMENDED THAT THE ORIGINAL .DAT FILE BE BACKED UP BEFORE USING THIS OPTION

Select **Import** to import data into DS-5. The data should either be a text file or a DS-5 file.

For text files the following format is required:

Dd/mm/yyyy

Hh:mm:ss

RECS_PER_SECOND (or SECONDS_PER_REC)

Nn

A

B

C

where the date is Dd/mm/yyyy.

the time is Hh:mm:ss (24 hour clock)

RECS_PER_SECOND means that the following number Nn represents the number of samples per second

SECONDS_PER_REC means that the following number Nn represents the number of seconds between each sample.

A, B, C represents the single data stream as one column


The data is picked up from the chosen file (date, time, duration and sample rate) and the entry boxes in the Import box allow the correction of the data entry for date, and time. The channel heading and units are entered e.g. Saturation and Pulse Rate.

Data may also be imported into a new file.

File | Import

Combine Files


Select **File | Combine Data files** to combine two or more DS-5 recordings and save as another file.



File | Combine Data Files

- **Files** - select the files by clicking on the first file and then, whilst holding the Ctrl or Shift key, click on the other files you wish to combine.
- **Directories** - shows the directory of files you wish to combine. The box allows navigation.
- **Files to combine** – Highlight and click the >> button to list files to be combined. Highlight and click the << button to remove files from the combined file list.
- **New file's time** shows start time, date and duration. Time between files is padded with zeros and the file duration increases accordingly.
- **New Name** - enter the title for the new file.
- **Combined Filename**: the filename for the new file is automatically generated by the program and is created using part of the filename of the first file selected. A different filename may be entered.

Analyse and Dump (Export) Patient information




File | Analyse and Dump patient info

Select Analyse and Dump patient info to export into a spreadsheet DIF file or an ASCII text file the basic patient and medical details along with analyses performed. The format is shown in Appendix A. The first line is the list of data exported, followed on the next lines by the data.

Combine Patient Information Files

Select Combine Patient Information Files to join together the spreadsheet DIF files or the ASCII text files which are exported by Analyse and Dump Patient Information




File | Combine Patient Information Files

- **Patient Information Files** - select the files you want to combine by clicking on the first file and then, whilst holding the Ctrl or Shift key, click on the other files.
- **Directories** - shows the directory of files you wish to combine. The box allows navigation.
- **Patient Information Files to combine** – click >> button to list files to be combined. Click << button to remove files from combined file list.
- **Combined Filename**: the filename for the new file is automatically generated by the program and is created using part of the filename of the first file selected. A different filename may be entered.

Group Auto Operation


It is possible to operate on one or more files at a time.



File | Group Auto Operation

- **Files** - select the files you want to operate on by clicking on the first file and then, whilst holding the Ctrl or Shift key, click on the other files.
- **Print Report** – analyse one or more files and print reports.
- **Analyse & Dump Info** – analyse one or more files and send joint patient information to the spreadsheet DIF file or the ASCII text file.
- **PDF Report** – analyse one or more files and create PDF reports.


Delete (a recording, analysis) and Delete Workpad

Click on the option **File | Delete** to delete an 'opened' recording or analysis or the toolbar icon **Delete** 

Recordings are sent to the Recycle bin from whence they may be recovered, but analyses (which may be recalculated) are deleted irrevocably.

Properties: (Information about the recording)


Select **File | Properties** to show information about the recording: includes the Title (name of the patient), start date and time of the recording, duration, source, averaging time and sampling rate of the different signals.



File | Properties

Window Menu

- Select the **Cascade** option to arrange opened documents (recordings or analyses) so that their title bars (containing names) are visible:
- Select the **Tile** option to tile vertically opened documents (recordings or analyses), so that they are visible, see below:




Window | Tile

- Select the **Arrange Icons** option to place/ arrange the minimised documents (as icons) into rows at the bottom of the current screen
- Select the **Close All** option to close the open documents (a save dialogue box will appear for each unsaved one).

List of opened documents

This section in the Windows menu displays the open documents in the order they were opened with the active document ticked.



Window choices

Standards

Safety and effectiveness

For use following consultation with a physician. The accuracy of the software is dependent upon the conformance and correct operation of the device as detailed in the operating manual.

Indications for use

DS-5 is a software programme which provides clinicians with the ability to collect, edit and analyse trend data made by devices.

It retrieves data recordings and provides general functions that are not intended to diagnose any specific disease.

DS-5 may be used to retrieve the stored results of the measurement by the device during patient monitoring.

User Responsibility

The user of the product has the responsibility for the correct installation and maintenance of the software, and the sole responsibility for any malfunction which arises from improper use, faulty maintenance, improper alteration, or damage by anyone other than KONICA MINOLTA SENSING, INC. authorised service personnel.

In particular, the manufacturers are not responsible for any malfunction caused by conflicting software.

Warnings and contra-indications

It is specifically required that the results be reviewed and interpreted by a clinician.

The accuracy of the recording and analyses may be affected by the position of the transducer and condition of the subject as well as the recording device.

Standards compliance

KONICA MINOLTA SENSING, INC. DS-5 meets EN 60601-1-1 and EN 60601-1-4. It is in risk category Class I according to Annex IX of the EC Council Directive 93/42/EEC of 14 June 1993.

Symbols

CE: The CE mark on this product denotes conformity with the EC Council Directive 93/42/EEC concerning medical devices.

Specifications

Measurement	None inherent, measurements are made by the monitor.
Range	Saturation 0-100%.
	Pulse Rate 0-250 bpm.

Other channels:

Accuracy	as data provided
Operating Conditions	not applicable
Storage Conditions	not applicable
Power	not applicable
Dimensions	typically 20Mb
Standards	EN 60601-1-1

Manufacturer's Address

KONICA MINOLTA SENSING, INC.
3-91, Daisennishimachi,
Sakai-ku, Sakai-shi,
Osaka,
590-8551,
Japan

Appendix A: Dump data

Spreadsheet DIF file

The first line is the list of data exported, followed on the next lines by the data, using a new line for each analysis e.g.

Filename	Name	Patient ID Number	Sex	Date of Birth	Age	Height	Weight	BMI
C:\DATA\Pulsox-300i\Pulsox-300i.dat	Sleep Apnoea, George, Mr	12345	M	01/12/1950	40	1.57	80	32.46

Excel spreadsheet (first few cells only) The analysis is in Row 2

ASCII text file

The first line is the list of data exported, with inverted commas delineating each, followed on the next lines by the data, using a new line for each analysis e.g.

"Filename" "Name" "Patient ID Number" "Sex" "Date of Birth" "Age" "Height" "Weight" "BMI" "Neck" "Epworth Score" "Address1" "Address2" "City" "County" "Post Code / ZIP" "Tel(Home)" "Tel(Work)" "Comments" "Health Authority" "Referring Physician" "GP" "GP Address1" "GP Address2" "GP City" "GP County" "GP Post Code / ZIP" "GP Tel" "Body No" "Date of Study" "Start Time" "Duration" "Label 1" "User data 1" "Label 2" "User data 2" "Label 3" "User data 3" "Label 4" "User data 4" "Label 5" "User data 5" "Label 6" "User data 6" "Label 7" "User data 7" "Label 8" "User data 8" "Label 9" "User data 9" "Label 10" "User data 10" "Label 11" "User data 11" "Label 12" "User data 12" "SpO2 Marker Mode" "SpO2 Analysis Interval" "Total time analysed for SpO2" "Total time analysed for apnoea index" "Events Rejected" "SpO2 cutoff for analysis" "SpO2 time below cutoff" "(a-b) - rise after drop" "c1 - Minimum time to drop in seconds" "c2 - Maximum time to drop in seconds" "d - Maximum time to rise in seconds" "Moving Average Time (sec)" "SpO2 Mean" "SpO2 Median" "SpO2 5th Percentile" "SpO2 95th Percentile" "Dip Level 1" "No. of Dips" "Dips/Hour" "Mean Dip nadir" "Dip Level 2" "No. of Dips" "Dips/Hour" "Mean Dip nadir" "Dip Level 3" "No. of Dips" "Dips/Hour" "Mean Dip nadir" "Min SpO2" "Time of Min SpO2" "12Sec Delta" "Delta Index" "Delta Index Sample Rate" "User value" "Time spent below user value" "Percentage of Time spent below user value" "Time spent SpO2 < 95%" "Percentage of Time spent SpO2 < 95%" "Time spent SpO2 < 90%" "Percentage of Time spent SpO2 < 90%" "Time spent SpO2 < 85%" "Percentage of Time spent SpO2 < 85%" "Time spent SpO2 < 80%" "Percentage of Time spent SpO2 < 80%" "Time spent SpO2 < 75%" "Percentage of Time spent SpO2 < 75%" "Time spent SpO2 < 70%" "Percentage of Time spent SpO2 < 70%" "Time spent SpO2 < 65%" "Percentage of Time spent SpO2 < 65%" "Time spent SpO2 < 60%" "Percentage of Time spent SpO2 < 60%" "Pulse Marker Mode" "Pulse Analysis Interval" "Total time analysed for Pulse Rate" "Total time analysed for Pulse Rises" "Events Rejected" "Pulse cutoff for analysis" "Pulse time below cutoff" "Pulse Mean" "Pulse Median" "Pulse Mode" "Pulse SD." "Pulse Definition [1]" "Pulse Rise" "Rises/Hour" "Pulse Definition [2]" "Pulse Rise" "Rises/Hour" "Pulse Definition [3]" "Pulse Rise" "Rises/Hour" "Max Pulse Peak" "Time of Max Pulse Peak" "C:\Documents and Settings\All Users\Application Data\KonicaMinolta\DS-5\Data\Pulsox-300i (Sleep Apnoea1)\Pulsox-300i (Sleep Apnoea1).dat" "Sleep Apnoea" "12345" "M" "" "040 01.57 80.00 32.46 "50.00" "10" "Address1" "Address2" "London" "UK" "123-45" "111-1111" "222-2222" "Comment #1 Comment" "Mr. H" "Dr. P" "Dr. G" "Address1" "Address2" "New York" "USA" "123-45" "333-3333" "0" "2005/05/24" "20:57:41" "08:00:00" "Label1" "Detail1" "Label2" "Detail2" "Label3" "" "Label4" "" "Label5" "" "Label6" "" "Label7" "" "Label8" "" "Label9" "" "Label10" "" "Label11" "" "Label12" "" "Full Study" "" "08:00:00" "07:29:17" "00:30:43" "00:30:37" "25" "00:00:06" "0" "4" "120" "20" "12" "91.62" "92.70" "81.50" "97.60" "4" "525" "65.64" "84.72" "3" "536" "67.01" "84.85" "2" "551" "68.89" "85.08" "62.90%" "0:36:52" "4.43" "5.89" "1" "< 88%" "01:33:11" "20.74%" "05:03:33" "67.56%" "02:23:39" "31.97%" "00:44:49" "9.98%" "00:17:45" "3.95%" "00:07:19" "1.63%" "00:01:40" "0.37%" "00:00:00" "0.00%" "00:00:00" "0.00%" "Full Study" "" "08:00:00" "07:29:17" "00:30:43"

Comparison between the patient information file of DS-5(Ver.1) and DS-5(Ver.2)

DS-5(Ver.1)	DS-5(Ver.2)
SpO2 Time Analysed	Total time analysed for SpO2
SpO2 Total time missed	Total time analysed for apnoea index
Pulse Time Analysed	Total time analysed for Pulse Rate
Pulse Total time missed	Total time analysed for Pulse Rises

Appendix B: KONICA MINOLTA Pulsox-3 series Oximeters: 3i/ 3iA /3Si/ 3Li Data Storage

The KONICA MINOLTA Pulsox-3 series oximeters store trend data as a new 'File' each time they are switched on.

The sampling rate for trend data is one sample every 5 seconds. Memory capacity is either 12 or 24 hours.

Oximeter time and date: These oximeters have internal clocks but these should be reset by the user from the computer if the battery power is removed for more than a short period. If the time and date have not been set, a default time and date (01/01/1997) will be set by the oximeter from when new batteries are inserted.

Clearing the oximeter memory should be routinely carried out after each recording is downloaded to the computer. This may be performed with the Clear Memory facility. If memory is not made available, the earliest data is overwritten.

Note: The Date/Time and File Clear and Output buttons on the IF-3 interface are only for use with the thermal printer which is not required if DS-5 software is used.

Communicating with the oximeter

Insert the oximeter in the Pulsox IF3 interface with nothing blocking the optical pathway between the oximeter and the sensor on the interface.

Connect the cable from the interface connector to the serial port of the computer, connect the IF3 interface to the AA-5P AC adapter, the adapter to mains power and slide the POWER switch to ON. Switch on the oximeter and wait until it has completed its startup test.

If the time and date have been set these values will appear in the startup sequence, followed by the F (File number) of the current recording.

Then

A flashing 'C' indicates that no probe is connected.

A flashing 'L' or 'P' indicates that the light or pulse signal is low.


If the low battery warning appears, new batteries should be inserted before continuing. The manual warns that replacement of the batteries must be completed within 5 minutes otherwise the measurement values and preset date/time will be lost.

Downloading Data

Pulsox-3i connected to


The drop down list shows the available 'COM' or 'serial' ports. Select the communication port to which the device is connected with the down arrow.

Choose Download, Set Clock or Clear Memory. (See sections following Downloading Data for Set Clock and Clear Memory)



Download Pulsox 3i dialogue box

On successful connection the LCD screen of the Pulsox-3i will display 'F out' and the information about recordings made is sent to the computer.



Download with several recordings on the Pulsox-3i recorder

Show All Files/ Hide Short Files Button

Short recordings of less than 300 seconds duration may be hidden from view or shown with this button.

To download recordings from the recorder to the computer

File Selection

Highlight one or more recordings. You can select two or more recordings by pressing <Shift> or <Ctrl> key simultaneously with the left mouse button.

Note: If you select two or more recordings DS-5 will automatically merge the recordings into one file with zero values added to fill the gaps between the recordings. These files may be very large (up to 14 days) so if you do not want to merge recordings they should be downloaded one at a time.

Location of downloaded file / Browse

'Location of downloaded file' shows where the selected recording will be saved on the computer. By default this is the program data directory.

To change where the file will be saved click the 'Browse' button and select a folder or create a new folder.

File Name

A unique File Name is automatically generated and the default location of the file is indicated.
A different filename may be entered and the location may be changed with the Browse button.

Surname / ID Number

You may enter the patient surname/ID to be saved in the file. If you leave the patient name field empty it will automatically be saved with a blank patient name, which could lead to problems in locating the file in the future.

Details

This gives access to a dialogue box to allow patient details to be entered. See 'Patient Details' section for the details.

Download Channels

Choose the channels to be downloaded by marking the check boxes.

Auto print report after downloading

You may select 'auto print report after downloading', to download, analyse and print the report in one operation.

Auto create PDF report after downloading

You may select 'auto create PDF report after downloading' to download, analyse and create the PDF report in one operation.

Auto create patient information file after downloading

You may select 'Auto create patient information file after downloading' to download, analyse and create the patient information file in one operation.

Report Setup

After downloading you may use the Analysis setup dialogue box to define the parameters used in the analysis. See 'Analysis: Run an Analysis and make a Report' section for the details.

Downloaded


When the OK button is clicked the highlighted recordings will be transferred to the computer. The data is saved by this action but it is not deleted from the recorder. The number of times a file has been downloaded is shown.

After the data has been transferred, the first few minutes will be shown on the 'Chart Screen'.

After downloading the data or selecting **Cancel** the Pulsox-3 series device is switched off. If the user wishes to carry out another task it must be switched back on.

Set Clock

Clicking on the **Set Clock** button will start the oximeter clock setting firmware. Upon successful connection a year (starting with '97 by default) will blink on the LCD screen of the oximeter and the 'set clock' dialogue box will be displayed on the computer.



Set Clock Dialogue Box

The text between the UP and DOWN buttons indicates the item of choice (initially this is 'Year'). These date/ time items are in sequence Year, Month, Day, Hour, Minute.

Changing a date/ time item

To increase or decrease an item, click the **UP** or **DOWN** buttons. The values will be shown on the oximeter LCD screen. Where there is more than one date/ time item shown on the LCD, eg month and day, the value of the chosen item will be shown flashing.

Moving to the next date/ time item

Click on the button labelled with the next available date/ time item to be moved to. The new values become valid after the Finish button is clicked. Previous time and dates set for existing stored recordings are not affected.

After clicking on Finish the Pulsox-3 series device is switched off. If the user wishes to carry out another task it must be switched back on.

Clear Memory

Clicking on the **Clear Memory** button and selecting the YES button in the warning dialogue box that follows will delete the oximeter memory.

The oximeter LCD will display '-- ---' for a second before the power is switched off.